
Nanoservice-enabled event-driven clinical decision support REST APIs


Hyponatremia Diagnosis Demo

For more information about Nanoservice, please visit https://csml.uconn.edu/

Step 1. Go to http://csmlab7.uconn.edu/swagger/




Step 2. [setValue] Click "[POST] /setValue"



Step 3. [setValue] Set the value of the variable "bloodSodium" to 98 by following the instruction.



- If the actor ID or variable does not exist, it will be automatically created.
- If "publish" is true, "actorId", "variable", and "value" are published to an Azure Event Grid topic (called "topic").



The code for publishing an event to **Azure Event Grid** https://azure.microsoft.com/en-us/services/event-grid/

```
public async Task PublishToAzureEventGridAsync(JObject data)
    string AzureEventGridTopicEndPoint = "https://topic.eastus-1.eventgrid.azure.net/api/events?api-version=2018-01-01";
    string AzureEventGridTopicAccessKey = "bHLip04YkH3Ysh0WvISAEUINVk3BWcPGTqGB6t/0iQw=";
   // [Warning] Be careful not to add "/" at the end of publisherBaseUri
    //string publisherBaseUri = "http://nanoservice3.eastus.cloudapp.azure.com";
    string publisherBaseUri = "http://csmlab7.uconn.edu";
    string uri = AzureEventGridTopicEndPoint;
    string topicSubject = (string)data.SelectToken("variable");
    data.Add("publisherBaseUri", publisherBaseUri);
   // Event data schema (Azure Event Grid)
    // https://docs.microsoft.com/en-us/azure/event-grid/post-to-custom-topic#event-data
                                                                            The source code is available at:
    dynamic requestBody = new ExpandoObject();
    requestBody.id = "notSet";
    requestBody.eventType = "notSet";
                                                                            https://github.com/yshin1209/Nanoservice
    requestBody.subject = topicSubject; // e.g., bloodSodium
    requestBody.eventTime = DateTime.Now;
    requestBody.data = data;
    requestBody.dataVersion = "v1";
    List<dynamic> requestBodyArray = new List<dynamic>();
    requestBodyArray.Add(requestBody);
    using (HttpClient client = new HttpClient())
       client.DefaultRequestHeaders.Add("aeg-sas-key", AzureEventGridTopicAccessKey);
       var request = new HttpRequestMessage(HttpMethod.Post, uri);
       string jsonRequestBody = JsonConvert.SerializeObject(requestBodyArray);
       request.Content = new StringContent(jsonRequestBody, Encoding.UTF8, "application/json");
       await client.SendAsync(request);
```

Azure Function https://azure.microsoft.com/en-us/services/functions/


```
public async static void Run(JObject eventGridEvent, TraceWriter log)
11 {
12
       log.Info(eventGridEvent.ToString(Formatting.Indented));
       var jsonEventGridEvent = eventGridEvent.ToString();
13
       JObject jobj = JObject.Parse (jsonEventGridEvent);
14
15
       double bloodSodiumValue = (double)jobj["data"]["value"];
       string actorId = (string)jobj["data"]["actorId"];
16
17
       if (bloodSodiumValue < 135) {decision = "true";}</pre>
18
       if (bloodSodiumValue >= 135) {decision = "false";}
19
20
       string publisherBaseUri = (string)jobj["data"]["publisherBaseUri"];
21
       string setValueUri = publisherBaseUri + "/setValue";
       var setRequest = new HttpRequestMessage(HttpMethod.Post, setValueUri);
       dynamic setRequestBody = new ExpandoObject();
24
               setRequestBody.variable = "hyponatremia";
26
               setRequestBody.value = decision;
27
       string jsonsetkequestbody - Jsonconvert.serializeouject(setRequestBody
28
       setRequest.Content = new StringContent(jsonSetRequestBody, Encoding.U)
       var secondResnonse = await client SendAsync(setRequest):
20
```

The Hyponatremia Diagnosis Function 1) decides whether there is hyponatremia or not, 2) adds a new variable "hyponatremia" to the patient actor in case the variable is not yet added and sets the value to true (hyponatremia) or false (not hyponatremia).

Step 4. [getValue] Get the value of the variable "hyponatremia" which is now "true" (because 98 < 135)

```
Example Value | Model


{
    "actorId": "patient032904475",
    "variable": "hyponatremia"
}
```




Step 5. [setValue] Set the value of the variable "bloodSodium" to 235.

```
{
    "actorId": "patient032904475",
        "variable": "bloodSodium",
        "value": 235,
        "publish": true
}
```

Step 6. [getValue] Get the value of the variable "hyponatremia" which is now "false" (because 235 > 135)



Step 7. [removeVariable] Remove the variable "bloodSodium".





For more information about Nanoservice, please visit https://csml.uconn.edu/